

BLEACH SACHETS

Bleach Powder Sachets



FEATURES

- 0.2% of available chlorine
- Easy to use
- Rapid disinfection
- Long shelf-life
- Convenient and stable source of chlorine

AREAS OF USE

- Institutional

DESCRIPTION

BLEACH SACHETS provide a safe and simple disinfecting agent for cleaning and sanitising the exterior and interior of syringes. BLEACH SACHETS feature a long shelf life, showing minimal chlorine loss after more than one year of storage. They come in individual 6g sachets, which provide 0.2% of available chlorine when diluted in 200mL of water.

DIRECTIONS FOR USE

Dilute one sachet in 200mL (3% w/v) of water to provide 0.2% (w/v) available chlorine. Further instructions are provided on the outside of each sachet.

Consult your Jasol representative for further information.

SAFETY INFORMATION

BLEACH SACHETS are classified as Hazardous according to the GHS. BLEACH SACHETS are classified as Dangerous Goods according to the ADG Code (7th Ed). BLEACH SACHETS are NOT a Scheduled Poison according to the SUSMP. Avoid contact with eyes and skin. Refer to SDS for further information.

FIRST AID INFORMATION

For advice in an emergency, contact a **Poisons Information Centre (Phone Australia 13 11 26 or New Zealand 0800 764 766) or a doctor at once.** If swallowed, do NOT induce vomiting. Wash out mouth thoroughly with water. If in eyes, hold eyelids apart and flood with water for at least 15 minutes. If on skin, remove contaminated clothing and wash affected area with running water. Seek immediate medical attention.

24 HOUR TOLL FREE EMERGENCY PHONE NO: 1800 629 953

PRODUCT ENQUIRIES

Customer Service: **1800 334 679**
jasolcsteam@gwf.com.au

For further information regarding product safety and handling, refer to the Safety Data Sheet and product label.

PHYSICAL & CHEMICAL PROPERTIES

FORM	Solid Granules
COLOUR	Pink
ODOUR	Chlorine
pH	~8.7 (1% Solution)
AVAILABLE CHLORINE	0.2% (w/v)
FLAMMABILITY	Not Combustible.

REFER OVERLEAF FOR ADDITIONAL INFORMATION



BLEACH SACHETS

Bleach Powder Sachets



INSTITUTIONAL / GOVT

EFFICACY

BLEACH SACHETS were evaluated in accordance with UNI EN 14476+A2:2019 at a concentration of 3% for 2 minutes at 20°C in hard water under clean conditions.

The product demonstrated virucidal activity against the following virus under the above conditions:

- BVDV strain (*Bovine viral diarrhea virus as surrogate for Hepatitis C virus (HCV)*)

