

LEMON BLEACH

Chlorinated Bleach



FEATURES

- Removes most organic based stains
- Suitable for bleaching mildew
- Rapid sanitising action
- Whitening effect
- Pleasant fragrance
- Hospital Grade Disinfectant ARTG 279787

AREAS OF USE

- Stain removal
- Laundry whitening
- Disinfection
- Deodorising

PRODUCT ENQUIRIES

Customer Service: **1800 334 679**
jasolcsteam@gwf.com.au

For further information regarding product safety and handling, refer to the Safety Data Sheet and product label.

DESCRIPTION

LEMON BLEACH is a sodium hypochlorite based bleach containing 4% w/w of available chlorine. Use on white and colourfast cottons to remove stains, such as fruit, vegetable, grass, ink, collar dirt, nappy stains, egg, coffee, tea and other. It will sanitise and deodorise baths, showers, basins, cloths, dishwashers, linoleum and plastic garbage bins.

DIRECTIONS FOR USE

It is recommended the spot test in an inconspicuous area be conducted before use.

Stain removal: Use at 50mL/L of hot water. Soak article in a solution for 5 minutes and rinse.

Whitening: Use at 5mL/L of wash water for bleaching in washing machines. Always dilute before use. Do not allow concentrated product to come in contact with clothes.

Sanitising/deodorising: Use at 250mL/L of water for wiping down work surfaces, grout cleaning, cleaning of garbage bins and similar works.

Consult your Jasol representative for further information.

SAFETY INFORMATION

LEMON BLEACH is classified as Hazardous according to the GHS.

LEMON BLEACH is NOT classified as Dangerous Goods according to the ADG Code (7th Ed).

LEMON BLEACH is a Scheduled Poison according to the SUSMP.

Avoid contact with eyes and skin. Refer to SDS for further information.

FIRST AID INFORMATION

For advice in an emergency, contact **Poisons Information Centre Australia (Phone 13 11 26)** or a **doctor at once**. If swallowed do NOT induce vomiting. Wash out mouth thoroughly with water. If in eyes, hold eyelids apart and flood with water for at least 15 minutes. If on skin, remove contaminated clothing and wash affected area with running water. Seek immediate medical attention.

24 HOUR TOLL FREE EMERGENCY PHONE NO: 1800 629 953

PHYSICAL & CHEMICAL PROPERTIES

FORM	Clear Liquid
COLOUR	Pale Yellow
ODOUR	Lemon Fragrance
pH	11.5-13.5
CHLORINE CONTENT	4% (Packing Date)
FLAMMABILITY	Not Flammable



REFER OVERLEAF FOR ADDITIONAL INFORMATION

LEMON BLEACH

Chlorinated Bleach



EFFICACY

LEMON BLEACH was evaluated by the AOAC Hard Surface Carrier Test 991.47, 48, 49 (modified) at 1:4 v/v concentration & 10 minutes contact time in standard hard water.

The product passed the test against all prescribed organisms:

- *Pseudomonas aeruginosa*
- *Staphylococcus aureus*
- *Salmonella choleraesuis*

LEMON BLEACH was evaluated in accordance with the TGA Test for Hospital grade Disinfectants (Option A) at 1:20 v/v concentration.

The product passed the test against the below prescribed organisms:

- *Escherichia coli*
- *Proteus vulgaris*
- *Pseudomonas aeruginosa*
- *Staphylococcus aureus*

