



SUPERSCENT T

Hospital Grade Disinfectant & Reodorant

FEATURES

- Highly concentrated
- Hard water compatible
- Safe on painted surfaces
- Phosphate free
- Meets TGA Hospital Grade Disinfectants requirements (AUST L 346550)

AREAS OF USE

- Hospitals
- Nursing homes
- Restaurants
- Cafeterias
- Hotels & motels
- Other public facilities

DESCRIPTION

SUPERSCENT T is a highly perfumed air freshener designed for odour control and disinfection in hospitals, restaurants, schools and general industries. It removes odour causing bacteria and helps achieve a high level of sanitation. SUPERSCENT T replaces unpleasant body, food and sewage odours with a particular fresh and long lasting Tropicana fragrance, creating a very pleasant atmosphere.

DIRECTIONS FOR USE

Air Freshener: Dilute 1:5 (cold water) into trigger bottle and spray.

Hospital Grade Disinfectant: Dilute 1:15. Pre-clean/wipe surface, apply solution, leave for 10 mins then rinse off. Store unused solution for up to 24 hours at room temperature.

Detergent Reodorant: Dilute 1:100 (hot water) for general purpose mopping floors.

Hard surface disinfectant only. Not to be used on skin.

To prevent microbial contamination, close containers immediately after use and do not add or insert any foreign materials or objects.

Consult your Jasol representative for further information.

SAFETY INFORMATION

SUPERSCENT T is classified as Hazardous according to the GHS.

SUPERSCENT T is NOT classified as Dangerous Goods according to the ADG Code (7th Ed).

SUPERSCENT T is NOT a Scheduled Poison according to the SUSMP.

Avoid contact with eyes and skin. Refer to SDS for further details.

FIRST AID INFORMATION

For advice in an emergency, contact a **Poisons Information Centre (Phone Australia 13 11 26 or New Zealand 0800 764 766) or a doctor at once**. If swallowed do NOT induce vomiting. Wash out mouth thoroughly with water and seek immediate medical attention. If in eyes, hold eyelids apart and flood with water for at least 15 minutes. If on skin, remove contaminated clothing and wash affected area with running water. Seek immediate medical attention.

PRODUCT ENQUIRIES

Customer Service: **1800 334 679**

jasolcsteam@gwf.com.au

For further information regarding product safety and handling, refer to the Safety Data Sheet and product label.

PHYSICAL & CHEMICAL PROPERTIES

FORM	Clear Liquid
COLOUR	Green
ODOUR	Citrus
pH	6.5-7.5
SPECIFIC GRAVITY	1.00
FLAMMABILITY	Not Flammable

REFER OVERLEAF FOR ADDITIONAL INFORMATION





SUPERSCENT T

Hospital Grade Disinfectant & Reodorant

EFFICACY

SUPERSCENT T was evaluated in accordance with the TGA Test for Hospital Grade Disinfectants (Option A) and passed the test against the following organisms when was applied at a dilution rate of 1:15:

- *Escherichia coli*
- *Proteus vulgaris*
- *Pseudomonas aeruginosa*
- *Staphylococcus aureus*

SUPERSCENT T was evaluated by the AOAC Hard Surface Carrier Test 991.47, 48, 49 (modified) at a dilution rate of 1:15 & 10 minutes contact time in standard hard water. The product passed the test against all prescribed organisms:

- *Pseudomonas aeruginosa*
- *Staphylococcus aureus*
- *Salmonella choleraesuis*

SUPERSCENT T was evaluated in accordance with TMCV 006 at a dilution rate of 1:30 & 10 minutes contact time at room temperature. The product passed the test against the following organism:

- *Murine Hepatitis Virus (MHV)-1 ATCC/VR-261 (TGA designated surrogate to justify COVID-19 efficacy claim).*

